

Identifyn Standards for Optical Microscopy Methods™Identifyn® Primary Antibody**Rabbit anti-Human BAP1 Polyclonal, IgG***Cat ID# - PA-RP 000012***Mouse anti-Human H2Ax Monoclonal, IgG1 - (1H1D9/B7)***Cat ID# - PA-MM 000002*Identifyn® Secondary Antibodies**Identifyn® SRM Alexa Fluor™ 594 Affinity Purified Goat Anti-Rabbit, IgG (H + L)***Cat ID# - SA 000046***Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Mouse, IgG (H + L)***Cat ID# - SA 000015*

BAP1 (BRCA1-associated protein 1) is a deubiquitinase with critical roles in DNA repair, chromatin dynamics, and cell cycle regulation. Its ability to interact with other proteins central to the DNA damage response is significant in maintaining genomic stability. BAP1 has several identified protein-protein interactions within the HR pathway, which are critical to DNA Damage Response (DDR) signaling.

BAP1 and BRCA1, BARD1, Rad51 and Others

BRCA1 is a tumor suppressor gene critical for homologous recombination (HR), a high-fidelity DNA repair pathway for double-strand breaks. BAP1's interaction with BRCA1 and other DNA repair proteins, such as BARD1 and Rad51, is linked to its deubiquitinating activity, suggesting that it regulates critical protein function during DNA repair. This relationship may influence BRCA1's stability, recruitment to DNA damage sites, or ability to mediate HR. This interaction underlines the importance of BAP1 in promoting accurate DNA repair and maintaining genomic stability.

BAP1 and ATM

ATM (ataxia-telangiectasia mutated) is a central kinase in the DNA damage response, primarily activated by DNA double-strand breaks. It phosphorylates various downstream targets to initiate DNA repair pathways, including homologous recombination and non-homologous end joining (NHEJ). While the direct interaction between BAP1 and ATM isn't as well-documented as other interactions, BAP1 may indirectly influence ATM's activity by affecting the stability or function of proteins involved in DDR. Deubiquitination (BAP1) and phosphorylation (ATM) may well affect how cells respond to DNA damage.

BAP1 and HERC2

HERC2 is an E3 ubiquitin ligase involved in regulating protein ubiquitination. It has a known role in genomic stability and is implicated in the DNA damage response. BAP1's interaction with HERC2 likely involves deubiquitination, suggesting that BAP1 could regulate HERC2's activity or modulate its effects on other proteins in the DNA repair pathways. This interaction could be critical in balancing ubiquitination and deubiquitination, impacting DNA repair processes and cell cycle control.

ASXL1/ASXL2

ASXL1 and ASXL2 (additional sex combs-like 1 and 2) are other key interactors of BAP1. This complex can influence chromatin remodeling, potentially affecting the accessibility of DNA repair proteins to sites of DNA damage.

Implications in DNA Repair

These interactions underscore BAP1's involvement in DNA repair pathways and the broader DNA damage response. By interacting with BRCA1, BAP1 may facilitate homologous recombination, while interactions with ATM and HERC2 highlight its role in broader DDR signaling. Critical interaction with ASXL1/ASXL2 highlights its importance in chromatin remodeling, potentially affecting the accessibility of DNA repair proteins to sites of DNA damage. BAP1's deubiquitinating activity can stabilize essential proteins, promote the proper assembly of repair complexes, and maintain chromatin structure, all of which are necessary for efficient DNA repair and genomic integrity.

Mutations in BAP1 that disrupt these interactions can increase genomic instability and are associated with various cancers, underscoring the clinical relevance of these protein interactions for understanding disease mechanisms and developing therapeutic strategies.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, DNA damage was induced by treatment with 2 μ M etoposide-spiked complete growth media for 24 hours, then fixed in -20 °C methanol and blocked with [Thermo Fisher Fish Serum Blocker \(cat#37527X3\)](#). The primary antibody, Rabbit anti-Human BAP1 Polyclonal, IgG, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3 times in 1X PBS. Identifyn® SRM Alexa Fluor™ 594 Affinity Purified Goat Anti-Rabbit, IgG (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5 times in 1X PBS.

The second primary antibody, Mouse anti-Human H2Ax Monoclonal, IgG1 - (1H1D9/B7), was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3 times in 1X PBS.

Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5 times in 1X PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI \(cat#P36962\)](#).

Microscope

Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63x/1.25 Oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:

Z Stack - (3D) - Panel A

Z-Stack = 48 Slices (4.7 μ m)

Channels = 3

Scaling (Per Pixel) = 0.143 μ m X 0.143 μ m X 0.100 μ m

Image size (Pixels) = 398 X 448

Image Size (Scaled) = 58.86 μ m X 64.00 μ m

Bit Depth = 14 Bit

Image Center Position = X: 50.76 mm, Y: 46.35 mm

Stage Position = X: 50.76 mm, Y: 46.35 mm
ROI Center Offset = X: 1.14 μ m, Y: 0.00 μ m

Z Stack - (3D) - Panel B
Z-Stack = 48 Slices (4.7 μ m)

Channels = 3
Scaling (Per Pixel) = 0.143 μ m X 0.143 μ m X 0.100 μ m
Image size (Pixels) = 178 X 139
Image Size (Scaled) = 25.43 μ m X 19.86 μ m
Bit Depth = 14 Bit
Image Center Position = X: 50.76 mm, Y: 46.35 mm
Stage Position = X: 50.76 mm, Y: 46.35 mm
ROI Center Offset = X: 1.14 μ m, Y: 0.00 μ m

594 -

Reflector = 45 Texas Red
Beam Splitter = 585
Filter Excitation Wavelength = 540 - 580
Filter Emission Wavelength = 593 - 668
Contrast Method = Fluorescence
Light Source = X-Cite Xylis Lamp Intensity @ 20.00%
Exposure Time = 200 ms
Excitation Wavelength = 590
Emission Wavelength = 618
Effective NA = 1.25
Imaging Device = AxioCam 807 Mono
Camera Adapter = 1x
Depth of Focus = 1.20 μ m
Binning Mode = 2,2

488 -

Reflector = 38 HE eGFP
Beam Splitter = 495
Filter Excitation Wavelength = 450 - 490
Filter Emission Wavelength = 500 - 550
Contrast Method = Fluorescence
Light Source = X-Cite Xylis Lamp Intensity @ 25.00%
Exposure Time = 250 ms
Excitation Wavelength = 493
Emission Wavelength = 517
Effective NA = 1.25
Imaging Device = AxioCam 807 Mono
Camera Adapter = 1x
Depth of Focus = 1.00 μ m
Binning Mode = 2,2

DAPI - Acquired but not shown in Data Analysis

Reflector = 49 DAPI
Beam Splitter = 395
Filter Excitation Wavelength = 335 - 383
Filter Emission Wavelength = 420 - 470

Contrast Method = Fluorescence
 Light Source = X-Cite Xylis Lamp Intensity @ 10.00%
 Exposure Time = 25 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.25
 Imaging Device = AxioCam 807 Mono
 Camera Adapter = 1x
 Depth of Focus = 0.90 μm
 Binning Mode = 2,2

Post Processing - SIM Processing

Apotome SIM Image Processing = Display Mode "Optical Sectioning" - enabled correction selected.
 Normalization, Phase Corrector, Fourier Filter, and Deconvolution **were NOT used.**

3D Image Processing - Panel A

3D Maximum Projection = Precise
 Appearance = Threshold 13% all channels, Ramp 0% all channels, Maximum 100% all channels
 Background = Black
 Light = Brightness 100%, Enabled Tone Mapping
 Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Clipping Image Process - Panel B

Clipping planes are introduced in Panel B to highlight the Rabbit anti-Human BAP1 Polyclonal, IgG, and the Mouse anti-Human H2Ax Monoclonal, IgG1 - (1H1D9/B7) within the nuclear compartment by adjusting or cutting away portions of the image using a triple-wedge cut, removing parts of the image in the X-Y, X-Z, and Y-Z planes.

X-Y = Activate was selected, textured opaque was chosen, and a white outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.
 Position of Clip = -9%

Clip X = 0°

Clip Z = 0°

X-Z = Activate was selected, textured opaque was chosen, and a green outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.
 Position of Clip = -24%

Clip Y = -45°

Clip Z = 0°

Y-Z = Activate was selected, textured opaque was chosen, and a red outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = -8%

Clip Y = -30°

Clip Z = 0°

Panel B shows a wedge-cut view displayed from the top with a slight backward rotation, highlighting the region of wedge removal and revealing the full depth of signal across the

z-stack. This orientation provides an alternative view of wedge geometry and confirms intranuclear localization of BAP1, as indicated by colocalization with H2AX.

Summary

Panel A shows a three-dimensional maximum-intensity projection of Rabbit anti-Human BAP1 Polyclonal, IgG detected with the Identifyn® Goat anti-Rabbit 594 secondary antibody (red), together with Mouse anti-Human H2AX Monoclonal, IgG1 – (1H1D9/B7) detected using the Identifyn® donkey anti-mouse 488 secondary antibody (green). Following etoposide-induced DNA damage, BAP1 demonstrates increased nuclear recruitment and overlap within H2AX-associated damage regions, consistent with activation of DNA double-strand break repair pathways and chromatin-associated DNA damage response signaling. The integrated volumetric projection demonstrates that BAP1 recruitment occurs broadly within the damaged nuclear compartment and supports its established role in homologous recombination-associated repair processes and chromatin remodeling following genotoxic stress.

Panel B shows a three-dimensional volume-intensity wedge projection with orthogonal planes (X–Y, white; X–Z, green; Y–Z, red), demonstrating volumetric nuclear colocalization of BAP1 and H2AX following etoposide-induced DNA damage. Orthogonal reconstruction confirms that the overlap persists throughout the nuclear volume rather than arising from a single focal plane. However, despite the improved optical sectioning and contrast enhancement provided by structured illumination microscopy (SIM), the reconstruction remains fundamentally diffraction-limited. As such, these data support nuclear recruitment and damage-associated colocalization of BAP1 in response to etoposide-induced stress, but do not provide sufficient spatial resolution to interpret nanoscale chromatin organization, direct molecular interaction, or discrete repair-complex architecture. The observed structures therefore represent biologically relevant damage-associated recruitment domains rather than definitive sub-diffraction repair assemblies, reinforcing the need for higher-resolution modalities when interpreting fine structural organization within DNA damage response microenvironments.

Image Correction

NONE

Image by E.F.Simon

Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.