

Identifyn Standards for Optical Microscopy Methods™Identifyn Primary Antibody**Rabbit anti-Human γ H2Ax (Serine 139) Recombinant Monoclonal, Mouse Fc, IgG (BLR014M)***Cat ID# - PA-RR 000005*Identifyn® Secondary Antibody**Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(ab')₂ (H + L)***Cat ID# - SA 000062*

Gamma H2AX (γ H2AX) is a variant of the histone H2A protein that plays a crucial role in the cellular response to DNA double-strand breaks (DSBs). Its primary function is as a marker of DNA damage, particularly in the context of double-strand breaks (DSBs). When DSBs occur in chromatin, γ H2AX is rapidly phosphorylated specifically at the serine 139 residue (referred to as Ser139 phosphorylation or γ -phosphorylation) near the site of the break. This phosphorylation event forms γ H2AX foci, visualized via immunofluorescence microscopy techniques as an aggregate of fluorescently labeled protein at the site of breaks.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2 μ M Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with [Thermo Fisher Fish Serum Blocker cat#37527X3](#). The primary antibody, Rabbit anti-Human γ H2Ax Recombinant Monoclonal, Mouse Fc, IgG (BLR014M), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(ab')₂ (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI cat#P36962](#).

Microscope

Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, AxioCam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:

Single Plane (2D)

Channels = 2

Scaling (Per Pixel) = 0.143 μ m X 0.143 μ m

Image size (Pixels) = 919 X 805

Image Size (Scaled) = 131.29 μ m X 115.00 μ m

Bit Depth = 14 Bit

Image Center Position = X: 62.87 μ m, Y: 42.95 μ m

Stage Position = X: 62.86 μ m, Y: 42.95 μ m

ROI Center Offset = X: 1.14 μ m, Y: 0.00 μ m

647 -

Reflector = 50 Cy 5

Beam Splitter = 660

Filter Excitation Wavelength = 625 - 655

Filter Emission Wavelength = 665 - 715
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @20.00%
 Exposure Time = 200 ms
 Excitation Wavelength = 653
 Emission Wavelength = 668
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 1.30 μm
 Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @6.00%
 Exposure Time = 12 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 0.90 μm
 Binning Mode = 2,2

Split View - Panel A

The first panel on the left features Rabbit anti-Human γH2Ax Recombinant Monoclonal, Mouse Fc, IgG (BLR014M) visualized with Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(a'')₂ (H + L), the middle panel features DAPI nuclear staining, and the last panel on the right features the merged image of all channels.

2.5D Image Processing

Zeiss Zen Software converts a two-dimensional image into a height map, known as a 2.5D Image. The most significant extension in the Z-direction represents the intensity values, creating a pseudo-three-dimensional image that displays the fluorescent localization or colocalization of the collected signal, along with its intensity in arbitrary units. Here we show the fluorescent intensity of the Human γH2Ax Recombinant Monoclonal, Mouse Fc, IgG2a primary antibody labeled with Identifyn® SRM Alexa Fluor™ 647 (red)

2.5D Display

Render Mode = Surface
 Grid distance = 1
 Palette and Show Faces at Side are both Selected

2.5D Display Options

Angle X = 63
 Angle Y = -0
 Z Scaling = 1.0
 Ambient = 50
 Specular = 20
 Shininess = 50
 Light Height = 2.0

Summary

Panel A - Split View illustrates the signal and localization of the Rabbit anti-Human γ H2Ax Recombinant Monoclonal, Mouse Fc, IgG2. The left section features the Rabbit anti-Human γ H2Ax Recombinant Monoclonal, Mouse Fc, IgG2a primary antibody labeled with Identifyn® SRM Alexa Fluor™ 647 (red). The middle image shows the DAPI (blue), and the right presents a merged image of the three channels. The merge demonstrates that gamma H2AX Mouse Fc is localized to the nucleus, as indicated by DAPI staining.

Panel B - Zeiss Zen Software converts a two-dimensional image into a height map, known as a 2.5D Image. The most significant extension in the Z-direction represents the intensity values, creating a pseudo-three-dimensional image that displays the fluorescent localization or colocalization of the collected signal, along with its intensity in arbitrary units. Here we show the fluorescent intensity of the Human γ H2Ax Recombinant Monoclonal, Mouse Fc, IgG2a primary antibody labeled with Identifyn® SRM Alexa Fluor™ 647 (red).

As widefield is the first in our series of the "step down" method, as expected and seen as highly consistent, the staining appears pan-nuclear with some concentrated aggregates commonly referred to as "foci." As we move through widefield SIM and Confocal Microscopy, much will remain the same; however, as we enter super-resolution microscopy, we will see that gamma H2AX (Mouse Fc) displays three-dimensional localization and structural details, as we have previously shown in gamma H2AX and its direct conjugates.

Image Correction

None

Image by E.F. Simon
 Image Analysis by B.T. Bennett
 Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.