

Identifyn Standards for Optical Microscopy Methods™Identifyn® Primary Antibody**Mouse anti-Human H2Ax Monoclonal, IgG1 - (1H1D9/B7)***Cat ID# - PA-MM 000002*Identifyn® Secondary Antibodies**Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(ab')₂ (H + L)***Cat ID# - SA 000062***Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Rabbit, F(ab')₂ (H + L)***Cat ID# - SA 000020*

H2AX, a variant of the histone protein H2A, has a primary function that relates to the organization and structure of chromatin, as well as the cellular response to DNA damage. Here are the essential functions and roles of H2AX:

Nucleosome Formation: H2AX, like other histone proteins, contributes to the formation of nucleosomes. A nucleosome is the fundamental unit of chromatin, consisting of DNA wrapped around a core of histone proteins (including two molecules each of H2A, H2B, H3, and H4). This organization plays a crucial role in compacting DNA within the cell nucleus and regulating gene expression by controlling access to DNA.

DNA Damage Response (DDR): One of H2AX's key functions, distinct from those of other histone proteins, is its function in the DNA damage response. When DNA suffers double-strand breaks (DSBs), H2AX is phosphorylated at a specific serine residue (serine 139), leading to the formation of γ-H2AX (gamma H2AX). This phosphorylation serves as an early signal of DNA damage.

Recruitment of DNA Repair Proteins: Phosphorylation of H2AX to form γ-H2AX creates a binding site for various proteins involved in DNA repair. It helps recruit proteins involved in the recognition, signaling, and repair of DNA damage, such as MDC1 (Mediator of DNA Damage Checkpoint 1), 53BP1 (p53-binding protein 1), BRCA1 (Breast Cancer 1), and others. These proteins work together to repair the damage and maintain genomic stability.

Chromatin Remodeling: Through its role in signaling and repair, H2AX plays a role in chromatin remodeling during the DNA damage response. The formation of γ-H2AX marks regions where chromatin structure is altered to allow access for repair proteins and facilitate the repair process.

H2AX's functions are fascinating studies in molecular biology. They involve a complex interplay of structural and signaling roles. Its structural role in chromatin and its specialized role in signaling and facilitating the DNA damage response, particularly in double-strand breaks, are of immense interest. Through its phosphorylation to form γ-H2AX, it acts as a key marker for DNA damage and plays an essential role in the recruitment and organization of proteins implicated in the DNA repair process.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with [Thermo Fisher, Fish Serum Blocker cat#37527X3](#). The primary antibody, Mouse anti-Human H2Ax Monoclonal, IgG1 -

(1H1D9/B7), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(ab')₂ (H + L) was incubated at room temperature for 1 hour at a dilution of 1:1000.

The second primary antibody, [Actin Recombinant Rabbit Monoclonal Antibody \(JF47-01\)](#), was then incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Rabbit, F(ab')₂ (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI cat#P36962](#).

Microscope

Zeiss Axio Imager.M1 Widefield, using a Plan-Apochromat 63x/1.25 Oil M27, AxioCam 807, and X-Cite Xylis Lamp, was used with the following parameters for each dye:

Single Plane (2D) - Z-plane 46 of 91 total - Panel A

Z Stack - (3D) - Panel D

Z-Stack = 91 Slices (9 µm)

Channels = 3

Scaling (Per Pixel) = 0.143 µm X 0.143 µm X 0.100 µm

Image Size (Pixels) = 1012 X 1081

Image Size (Scaled) = 144.57 µm X 154.43 µm

Bit Depth = 14 Bit

Image Center Position = X: 49.25 mm, Y: 50.66 mm

Stage Position = X: 49.25 mm, Y: 50.66 mm

ROI Center Offset = X: 1.14 µm, Y: 0.00 µm

647 -

Reflector = 50 Cy 5

Beam Splitter = 660

Filter Excitation Wavelength = 625 - 655

Filter Emission Wavelength = 665 - 715

Contrast Method = Fluorescence

Light source = X-Cite Xylis Lamp Intensity @8.00%

Exposure Time = 80 ms

Excitation Wavelength = 653

Emission Wavelength = 668

Effective NA = 1.25

Imaging Device = 807 mono

Camera Adapter = 1X

Depth of Focus = 1.30 µm

Binning Mode = 2,2

488 -

Reflector = 38 HE eGFP

Beam Splitter = 495
 Filter Excitation Wavelength = 450 - 490
 Filter Emission Wavelength = 500 - 550
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @10.00%
 Exposure Time = 150 ms
 Excitation Wavelength = 493
 Emission Wavelength = 517
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 1.00 μm
 Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
 Beam Splitter = 395
 Filter Excitation Wavelength = 335 - 383
 Filter Emission Wavelength = 420 - 470
 Contrast Method = Fluorescence
 Light source = X-Cite Xylis Lamp Intensity @6.00%
 Exposure Time = 18 ms
 Excitation Wavelength = 353
 Emission Wavelength = 465
 Effective NA = 1.25
 Imaging Device = 807 mono
 Camera Adapter = 1X
 Depth of Focus = 0.90 μm
 Binning Mode = 2,2

Post Processing - SIM Processing

Apotome SIM Image Processing = Display Mod" "Optical Sectioning" - enabled correction selected.

Normalization, Phase Corrector, Fourier Filter, and Deconvolution **were NOT used.**

Split View - Panel A

The first panel on the left features Mouse anti-Human H2Ax Monoclonal, IgG1 - (1H1D9/B7) visualized with Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, F(ab')₂ (H + L), the second panel shows Actin visualized with Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Rabbit, F(ab')₂ (H + L), the third panel features DAPI nuclear staining. The fourth panel features the merged image of all channels. This split View was taken at Z-Position 46 of 91.

Profile View Display - Panel C

Profile Definition = Stroke Thickness 1, Normal View, Grid Distance 1 Profile View = X and Y axes were set to Auto - X (Min 0.0 and Max 105.1), Y (Min 0 and Max 668)

3D Image Processing - Panel D

3D Maximum Projection = Precise

Appearance = Threshold 13%, Ramp 0%, Maximum 100%

Background = Black

Light = Brightness 100%, Tone Mapping Enabled

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Clipping Image Process = Clipping planes are introduced to pronounce regions of special interest.

X-Y = Active was selected, textured opaque was chosen to remove a horizontal portion of the upper Z-Stack of the cell from the image, and a red outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = 5%

Clip X = 0°

Clip Y = 0°

Image Correction

None

Image by E.F. Simon

Image Analysis by B.T. Bennett

Data Presentation by P.D. Amistoso

© Copyright 2026 GMD12 LLC.