

Identifyn Standards for Optical Microscopy Methods™

Identifyn® Primary Antibody

Mouse anti-Human, ATRIP Monoclonal, IgG2b - (4H6H4/G3)

Cat ID# - PA-MM 000003

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, IgG (H + L), (Subclasses 1+2a+2b+3), Fcy Fragment Specific

Cat ID# - SA 000056

ATRIP (ATR Interacting Protein) is a key partner and co-activator of ATR (Ataxia Telangiectasia and Rad3-related), a crucial kinase in the DNA damage response, particularly in response to replication stress and single-strand breaks (SSBs). ATRIP plays a pivotal role in recognizing and binding to single-stranded DNA (ssDNA), facilitating ATR activation at DNA damage sites.

Overview of ATRIP

ATRIP is the regulatory subunit that associates with ATR and directs it to sites of DNA damage. It contains a domain that allows it to bind to Replication Protein A (RPA), a complex that binds and stabilizes single-stranded DNA. This interaction with RPA-coated ssDNA is crucial for ATRIP's ability to guide ATR to regions where the DNA damage response is needed.

Role of ATRIP in Single-Strand Break Repair

ATRIP plays a central role in the activation of ATR, which is critical to the cellular response to single-strand breaks and replication stress:

Recognition of Single-Stranded DNA: ATRIP, through its interaction with RPA, identifies regions of ssDNA. These regions are commonly found at stalled replication forks, resected DNA ends, or sites of single-strand breaks.

Recruitment of ATR: Once ATRIP binds to RPA-coated ssDNA, ATR is recruited to the site. This recruitment is crucial for ATR activation, leading to a cascade of events that coordinate the DNA damage response and facilitate DNA repair.

Activation of Checkpoints: By facilitating ATR activation, ATRIP indirectly contributes to checkpoint signaling. This leads to CHK1 phosphorylation, which plays a role in cell cycle arrest, providing time for the repair of single-strand breaks before the cell cycle progresses.

ATRIP's Broader Role in the DNA Damage Response

Beyond its specific role in single-strand break repair, ATRIP's interaction with ATR has broader implications for the DNA damage response:

Coordination of DNA Repair: The ATR/ATRIP complex orchestrates the recruitment of DNA repair proteins. This coordination ensures that repair processes like base excision repair (BER) and homologous recombination can proceed efficiently at sites of single-strand breaks.

Response to Replication Stress: ATRIP is crucial in responding to replication stress, which can lead to single-strand breaks or stalled replication forks. By facilitating ATR activation, ATRIP helps stabilize replication forks and prevent more severe DNA damage, such as double-strand breaks.

ATRIP is a critical regulatory protein that guides ATR to sites of single-strand breaks and other forms of replication stress. By recognizing RPA-coated ssDNA and recruiting ATR, ATRIP plays a pivotal role in activating the DNA damage response, leading to checkpoint activation and coordination of DNA repair processes. Its role in stabilizing replication forks and promoting genomic stability underscores its importance in maintaining genome integrity.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2µM Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with [Thermo Fisher Fish Serum Blocker cat#37527X3](#). The primary antibody, Mouse anti-Human, ATRIP Monoclonal, IgG2b - (4H6H4/G3), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Mouse, IgG (H + L), (Subclasses 1+2a+2b+3), Fcy Fragment Specific was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI cat#P36962](#).

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1.25X Tubelens, was used with the following parameters for each dye:

Single Plane (2D) - Z-plane 63 of 139 total - Panels A

Channels = 2
Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm
Image Size (Pixels) = 1514 X x 1368
Image Size (Scaled) = 31.96 µm X 28.88 µm
Bit Depth = 16 Bit
Image Center Position = X: 62.03 mm, Y: 38.74 mm
Stage Position = X: 62.03 mm, Y: 38.74 mm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

Z Stack - (3D) - Panel B

Z-Stack = 139 Slices (4.14 µm)

Channels = 2
Scaling (Per Pixel) = 21.11 nm X 21.11 nm X 30.00 nm
Image Size (Pixels) = 234 X x 228
Image Size (Scaled) = 4.94 µm X 4.81 µm
Bit Depth = 16 Bit
Image Center Position = X: 62.03 mm, Y: 38.74 mm
Stage Position = X: 62.03 mm, Y: 38.74 mm
ROI Center Offset = X: 0.00 µm, Y: 0.00 µm

647 -

Reflectors = SR HE LBF 405/488/561/642
Beam Splitter = 405, 488, 561, 642
Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method = Fluorescence
Illumination Wavelength = 640
Channel Name = T1-Axiocam 820-SR
Excitation Wavelength = 640
Emission Wavelength = 729
Effective NA = 1.4
Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Imaging Device = Axiocam 820
Camera Adapter = 1X
Exposure Time = 120 ms
Depth of Focus = 1.13 μ m
Binning Mode = 2,2
Laser Wavelength = 640 nm @1.00%
Frame Time = 1.59 s
Scan Direction = Unidirectional
Grating = 36.5 μ m grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
Beam Splitter = 405, 488, 561, 642
Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method = Fluorescence
Illumination Wavelength = 405
Channel Name = T3-Axiocam 820 #2-SR
Excitation Wavelength = 405
Emission Wavelength = 445
Effective NA = 1.4
Detection Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Imaging Device = Axiocam 820 #2
Camera Adapter = 1X
Exposure Time = 50 ms
Depth of Focus = 0.69 μ m
Binning Mode = 2,2
Laser Wavelength = 405 nm @2.00%
Frame Time = 676.00 ms
Scan Direction = Unidirectional
Grating = 36.5 μ m grating

Post Processing - SIM2 Processing

Processing Dimension - 3D
Result Sampling = 4
Process Sampling = 4

Channels = 2
Algorithm = SIM2
Input SNR = Low
Regularization Weight = 0.0000
Iterations = 3
Filter = No Filter
Fitting = Fast Fit
Normalization = Auto

3D Image Processing - Panel B

3D Volume Projection = Precise
Appearance = Threshold 2% all channels, Ramp 98% all channels, Maximum 100% all channels
Background = Black
Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection = View Angle 35°, Scale Z 30%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Summary

Using lattice SIM² super-resolution imaging, features that appear as single foci in Confocal and Airyscan images are resolved into distinct structures. However, only limited additional structural detail is gained when progressing from SIM to SIM². In contrast, single-molecule STORM imaging is required to more accurately resolve fine structural features of the ATRIP signal. Notably, while STORM provides enhanced spatial detail, additional organization likely remains below the resolution limit of this modality.

Image Corrections

DAPI was shown in white in panel B, rather than the default blue, to enhance visual contrast.

Image by P. Raut
Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.