

Identifyfyn Standards for Optical Microscopy Methods™Identifyfyn Primary Antibody**Rabbit anti-Human RAD50 Polyclonal, IgG***Cat ID# - PA-RP 000007*Identifyfyn® Secondary Antibodies**Identifyfyn®SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Rabbit, IgG (H + L)***Cat ID# - SA 000057***Identifyfyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Mouse, IgG (H + L)***Cat ID# - SA 000015*

Rad50 recognizes and processes double-strand breaks during DNA repair. It is a component of the MRN complex (Mre11-Rad50-Nbs1), which is crucial for identifying and processing DNA double-strand breaks. Rad50 also provides tethering and bridging of DNA ends. In its homodimeric state, Rad50 contains an ATPase domain at its ends, facilitating the alignment of broken DNA strands. Rad50 is also associated with telomere maintenance. Dysregulation or mutations in Rad50 or the MRN complex can lead to genomic instability, increased cancer susceptibility, and other genetic disorders.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency, then fixed in -20°C methanol and blocked with [Thermo Fisher Fish Serum Blocker \(cat#37527X3\)](#). The primary antibody, Rabbit anti-Human RAD50 Polyclonal, IgG, was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyfyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Rabbit, IgG (H + L), was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 3 times in PBS.

The second primary antibody, [Alpha Tubulin Monoclonal Antibody \(DM1A\)](#), was incubated for 1 hour at room temperature at a 1:250 dilution. Cells were washed 3X in PBS. Identifyfyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Mouse, IgG (H + L) was incubated for 1 hour at room temperature at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI, cat#P36962](#).

Microscope

Zeiss Axio Observer.Z1/7 Apotome SIM, using an EC Plan-Neofluar 63X/1.25 oil M27, Axiocam 807, and X-Cite Xylis Lamp was used with the following parameters for each dye:

Z Stack - (3D) - Panel A**Z-Stack = 89 Slices (8.8 µm)**

Channels = 3

Scaling (Per Pixel) = 0.143 µm X 0.143 µm X 0.100 µm

Image size (Pixels) = 1002 X 887

Image Size (Scaled) = 143.14 μm X 126.71 μm
Bit Depth = 14 Bit
Image Center Position = X: 87.68 μm , Y: 53.20 μm
Stage Position = X: 87.68 μm , Y: 53.20 μm
ROI Center Offset = X: 1.14 μm , Y: 0.00 μm

Single Plane (2D) - Panel B, Z-Plane 43 of 89, (Dimension Shown in Panel A)

Z Stack - (3D) - Panel C

Z-Stack = 89 Slices (8.8 μm)

Channels = 3
Scaling (Per Pixel) = 0.143 μm X 0.143 μm X 0.100 μm
Image size (Pixels) = 365 X 375
Image Size (Scaled) = 52.14 μm X 53.57 μm
Bit Depth = 14 Bit
Image Center Position = X: 87.68 μm , Y: 53.20 μm
Stage Position = X: 87.68 μm , Y: 53.20 μm
ROI Center Offset = X: 1.14 μm , Y: 0.00 μm

647 -

Reflector = 50 Cy 5
Beam Splitter = 660
Filter Excitation Wavelength = 625 - 655
Filter Emission Wavelength = 665 - 715
Contrast Method = Fluorescence
Light source = X-Cite Xylis Lamp Intensity @15.00%
Exposure Time = 100 ms
Excitation Wavelength = 653
Emission Wavelength = 668
Effective NA = 1.25
Imaging Device = Axiocam 807 Mono
Camera Adapter = 1X
Depth of Focus = 1.30 μm
Binning Mode = 2,2

488 -

Reflector = 38 eGFP
Beam Splitter = 495
Filter Excitation Wavelength = 450 - 490
Filter Emission Wavelength = 500 - 550
Contrast Method = Fluorescence
Light source = X-Cite Xylis Lamp Intensity @20.00%
Exposure Time = 200 ms
Excitation Wavelength = 493
Emission Wavelength = 517
Effective NA = 1.25
Imaging Device = Axiocam 807 Mono
Camera Adapter = 1X

Depth of Focus = 1.00 μm
Binning Mode = 2,2

DAPI -

Reflector = 49 DAPI
Beam Splitter = 395
Filter Excitation Wavelength = 335 - 383
Filter Emission Wavelength = 420 - 470
Contrast Method = Fluorescence
Light source = X-Cite Xylis Lamp Intensity @10.00%
Exposure Time = 25 ms
Excitation Wavelength = 353
Emission Wavelength = 465
Effective NA = 1.25
Imaging Device = AxioCam 807 Mono
Camera Adapter = 1X
Depth of Focus = 0.90 μm
Binning Mode = 2,2

Post Processing - SIM Processing

Apotome SIM Image Processing = Display Mode "Optical Sectioning" - enabled correction selected.

Normalization, Phase Corrector, Fourier Filter, and Deconvolution **were NOT used.**

3D Image Processing - Panel A

3D Maximum Projection = Precise
Appearance = Threshold 13% all channels, Ramp 0% all channels, Maximum 100% all channels
Background = Black
Light = Brightness 100%, Enabled Tone Mapping
Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

3D Image Processing - Panel C

3D Volume Projection = Precise
Appearance = Threshold 13% all channels, Ramp 50% all channels, Maximum 100% all channels
Background = Black
Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Directional Lighting, Enabled Tone Mapping
Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Clipping Image Process - Panel C

Clipping planes are introduced in Panel B to highlight the RAD50 Rabbit Polyclonal within the nuclear compartment by adjusting or cutting away portions of the image using a triple-wedge cut, removing parts of the image in the X-Y, X-Z, and Y-Z planes.

X-Y = Activate was selected, textured opaque was chosen, and a green outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = -3%

Clip X = 0°

Clip Z = 0°

X-Z = Activate was selected, textured opaque was chosen, and a red outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = -33%

Clip Y = 45°

Clip Z = 0°

Y-Z = Activate was selected, textured opaque was chosen, and a blue outline of the clipped perimeter was shown; Clip Back was selected, and Clip Front was NOT selected.

Position of Clip = -14%

Clip Y = -45°

Clip Z = 0°

Summary - RAD50 Recruitment Assessed by Widefield SIM (Apotome)

Three-dimensional maximum-intensity projection showing RAD50 detected with an Identifyn® polyclonal primary antibody and visualized using Identifyn® SRM Alexa Fluor™ 647 goat anti-rabbit secondary antibody (pseudo-colored white). α -Tubulin is visualized with an Identifyn® SRM Alexa Fluor™ 488 donkey anti-mouse secondary antibody (green) and serves as a cytoplasmic exclusionary control, while DAPI marks nuclear DNA (blue).

ROI Panel C corresponds to Z-plane 43 of 89. A three-dimensional volume reconstruction (wedge view) illustrates orthogonal planes (X-Y, X-Z, Y-Z) of RAD50 following etoposide-induced DNA damage. Across all axes, the RAD50 signal is confined to the nuclear compartment, with no detectable cytoplasmic or nucleolar localization. Imaging was performed using low-intensity LED excitation and short exposure times to preserve structural fidelity and avoid acquisition-driven signal amplification.

Acquired using widefield structured illumination microscopy (Apotome), this dataset demonstrates that three-dimensional imaging at SIM resolution enables visualization of RAD50 aggregates that are not discernible in conventional widefield imaging. While these aggregates do not resolve the substructures observed at higher-resolution modalities downstream in the step-down workflow (SIM² and STORM), the signal provides a clear go / no-go readout confirming that etoposide induces a measurable RAD50 DNA damage response.

Image Correction

Identifyn® SRM Alexa Fluor™ 647 Affinity Purified Goat Anti-Rabbit, IgG (H + L) was pseudo-colored in Zen Software from red, the default color, to white, to provide visual contrast in Panels B and C.

Image by E.F.Simon

Image Analysis and Data Presentation by B.T. Bennett

© Copyright 2026 GMD12 LLC.

