

Identifyn Standards for Optical Microscopy Methods™

Identifyn Primary Antibody

Rabbit anti-Human γ H2Ax (Serine 139) Recombinant Monoclonal, Human Fc, IgG (BLR014M)

Cat ID# - PA-RR-000004

Identifyn® Secondary Antibodies

Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Human IgG (H + L)

Cat ID# - SA 000016

Gamma H2AX (γ H2AX) is a variant of the histone H2A protein that plays a crucial role in the cellular response to DNA double-strand breaks (DSBs). Its primary function is as a marker of DNA damage, particularly in the context of double-strand breaks (DSBs). When DSBs occur in chromatin, γ H2AX is rapidly phosphorylated specifically at the serine 139 residue (referred to as Ser139 phosphorylation or γ -phosphorylation) near the site of the break. This phosphorylation event forms γ H2AX foci, visualized via immunofluorescence microscopy techniques as an aggregate of fluorescently labeled protein at the site of breaks.

Method

HeLa cells were grown on acid-etched #1.5H cover glass to ~50% confluency. DNA damage was induced via treatment with 2 μ M Etoposide-spiked complete growth media for 24 hours, followed by -20°C methanol fixation and blocked with [Thermo Fisher Fish Serum Blocker cat#37527X3](#). The primary antibody, Rabbit anti-Human γ H2Ax Recombinant Monoclonal, Human Fc, IgG, was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS. Identifyn® 488 secondary antibody was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI cat#P36962](#).

Microscope

Zeiss LSM 900 AiryScan, Axio Observer.Z1/7, using with a Plan-Apochromat 63X/1.40 oil DIC M27 was used with the following parameters for each dye:

Z Stack - (3D) - Panel A

Z-Stack = 39 Slices (1.52 μ m)

Channels = 2

Scaling (Per Pixel) = 0.071 μ m X 0.071 μ m X 0.071 μ m

Image Size (Pixels) = 730 X 586

Image Size (Scaled) = 51.54 μ m X 41.37 μ m

Bit Depth = 16 Bit

Image Center Position = X: 77.00mm, Y: 53.10mm

Stage Position = X: 77.00mm, Y: 53.10mm

ROI Center Offset = X: 0.00 μ m, Y: 0.00 μ m

488 -

Reflector = None

Contrast Method = Fluorescence

Pinhole = 1.00 AU / 45 μ m

LMS MBS = MBS 405 + 488 + 561 + 640(T10/R90)

Laser Wavelength = 488nm @1.0%
Laser Blanking = Enabled
Scan Mode = Frame
Scan Zoom = 1.4
Rotation = 0°
Sampling = 1.00
Pixel Time = 0.51µs
Frame Time = 10.15s
LSM Scan Speed = 9
Scan Direction = Bidirectional
Line Step = 1
Averaging = 16
Excitation Wavelength = 493
Emission Wavelength = 517
Effective NA = 1.4
Detection Wavelength = 484-590
Imaging Device = Airyscan
Detector Type = GaAsP-PMT
Detector Gain = 650 V
Detector Offset = 0
Detector Digital Gain = 1.0

DAPI -

Reflector = None
Contrast Method = Fluorescence
Pinhole = 1.00 AU / 37µm
LMS MBS = MBS 405 + 488 + 561 + 640(T10/R90)
Laser Wavelength = 405nm @0.3%
Laser Blanking = Enabled
Scan Mode = Frame
Scan Zoom = 1.4
Rotation = 0°
Sampling = 1.00
Pixel Time = 0.51µs
Frame Time = 10.15s
LSM Scan Speed = 9
Scan Direction = Bidirectional
Line Step = 1
Averaging = 16
Excitation Wavelength = 353
Emission Wavelength = 465
Effective NA = 1.4
Detection Wavelength = 410-470
Imaging Device = GaAsP-Pmt1
Detector Type = GaAsP-PMT
Detector Gain = 650 V
Detector Offset = 0
Detector Digital Gain = 1.0

3D Image Processing - Panel A

3D Maximum Projection = Precise

Appearance = Threshold 1% (488), Threshold 1% (405), Ramp 0.0% all channels,

Maximum 100% all channels

Background = Black

Light = Brightness 100%, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Colocalization View - Panels B and C

Horizontal Axis - 488, Threshold 1475, Range Auto

Vertical Axis - DAPI, Threshold 384, Range Auto

Profile View Display

Profile Definition = Stroke Thickness 1, Normal View, Grid Distance 1

Profile View = X and Y axes were set to Auto - X (Min 0.00 and Max 13.4), Y (Min 4854 and Max 53425)

Summary

Panel A presents a three-dimensional view of the complete image acquisition, displayed as a maximum projection. An X-Y cut plane is positioned at Z plane 19 of 39 to emphasize the center of the nucleus and the gamma-H2AX signal within the nuclear compartment. Panels B and C demonstrate that the gamma-H2AX signal is localized within the nuclear compartment, as indicated in the colocalization plot. Panels D and E display a line profile across several foci, highlighting the peak fluorescent intensity and the distance between the foci signals.

It is important to note that a single Foci can be misleading at diffraction-limited confocal resolutions since it measures approximately 1 micron in width. In contrast, more advanced super-resolution techniques, such as STORM, can achieve a resolution of foci down to approximately 40 nm. This suggests our observed signal may represent more complex structures beyond our current resolution.

Image Corrections -

None

Image by P. Raut

Analysis and Presentation by B. T. Bennett

Copyright 2026 GMD12 LLC.