

Identifyn Standards for Optical Microscopy Methods™**Identifyn® SRM Alexa Fluor™ 488 Affinity Purified Donkey Anti-Mouse IgG (H + L)***Cat ID# - SA 000015*

Hela Cells were grown on acid-etched #1.5H cover glass to ~50% confluency, followed by -20°C methanol fixation and blocked with [Thermo Fisher Fish Serum Blocker cat#37527X3](#). The primary antibody, [Alpha Tubulin Monoclonal Antibody \(DM1A\)](#), was incubated for 1 hour at room temperature at a dilution of 1:250. Cells were washed 3X in PBS, and Identifyn® 488 secondary antibody was incubated at room temperature for 1 hour at a dilution of 1:1000. Cells were washed 5X in PBS and mounted using [ProLong™ Diamond Antifade Mountant with DAPI cat#P36962](#).

Microscope

Zeiss Elyra 5 Lattice SIM, AxioObserver.Z1/ 7 using Plan-Apochromat 63X/1.40 oil DIC M27, 1X Tubelens, was used with the following parameters for each dye:

Z Stack - (3D)**Z-Stack = 47 Slices (2.3µm)**

Channels = 2

Scaling (Per Pixel) = 16.89nm X 16.89nm X 50.00nm

Image Size (Pixels) = 1385 X 1359

Image Size (Scaled) = 23.39µm X 22.95µm

Bit Depth = 16 Bit

Image Center Position = X: 58.06mm, Y: 39.13mm

Stage Position = X: 58.06mm, Y: 39.13mm

ROI Center Offset = X: 0.00µm, Y: 0.00µm

488 -

Reflectors = SR HE LBF 405/488/561/642

Beam Splitter = 405, 488, 561, 642

Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642

Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800

Contrast Method = Fluorescence

Illumination Wavelength = 488

Channel Name = T2-Axiocam 820 -SR

Excitation Wavelength = 488

Emission Wavelength = 525

Effective NA = 1.4

Detection Wavelength = 419-470, 503-546, 576-617, 657-800

Imaging Device = Axiocam 820

Camera Adapter = 1X

Exposure Time = 100ms

Depth of Focus = 0.81µm

Binning Mode = 2,2

Laser Wavelength = 488nm @1.00%

Frame Time = 1.33s

Scan Direction = Unidirectional
Grating 32µm grating

DAPI -

Reflectors = SR HE LBF 405/488/561/642
Beam Splitter = 405, 488, 561, 642
Filter Excitation Wavelength = 200 - 405, 200 - 488, 200 - 561, 200 - 642
Filter Emission Wavelength = 419 - 470, 503 - 546, 576 - 617, 657 - 800
Contrast Method = Fluorescence
Illumination Wavelength = 405
Channel Name = T2-Axiocam 820 #2-SR
Excitation Wavelength = 405
Emission Wavelength = 445
Effective NA = 1.4
Detection Wavelength = 419-470, 503 - 546, 576 - 617, 657-800
Imaging Device = Axiocam 820#2
Camera Adapter = 1X
Exposure Time = 50ms
Depth of Focus = 0.69µm
Binning Mode = 2,2
Laser Wavelength = 405nm @7.0%
Frame Time = 676.00ms
Scan Direction = Unidirectional
Grating 36.5µm grating

Post Processing - SIM2 Processing

Processing Dimension - 3D
Result Sampling = 4
Process Sampling = 4
Channels = 1
Algorithm = SIM2
SIM Settings = Standard
Auto Sharpness = On
Sharpness = 10.3380
Input SNR = N/A
Regularization Weight = N/A
Iterations = N/A
Filter = N/A
Fitting = N/A
Normalization = Automatic
Experimental PSF = N/A

3D Image Processing

3D Volume Projection = Precise
Appearance = Threshold 2% (488), Threshold 2% (DAPI), Ramp 98% all channels,
Maximum 100% all channels
Background = Black
Light = Brightness 100%, Azimuth 80°, Elongation -130°, Enabled Tone Mapping

Projection = View Angle 35°, Scale Z 100%, Reverse Z-Axis NOT Selected, Stereo Anaglyph NOT Selected, Camera Separation 5% (deselected), Parallax Shift 0% (deselected)

Movie Series Settings

Render Series = Turn Around Y

Frames per Second = 120

Partial Panorama = Start Angle 0.0°, Stop Angle 90.0°

Image Corrections

None

Image by P. Raut

Copyright 2026 GMD12 LLC.

